

CURONIX

Continuing Your Pain Relief Journey

Post-Procedure Guide



We are excited to be a part of your pain relief journey with the Curonix Freedom® System. We are here to support you in connection with your device related issues.

We have provided answers to some of the frequently asked questions related to your new device for reference and informational purposes. Contact your local Curonix representative or Curonix Support at (800) 965-5134 for any device related questions or issues. Any medical questions should be directed to your Health Care Provider.





Table of Contents

Quick Reference User Guide

Transmitter Assembly	4-5
PNS Wearable	6
Torso Wearable	7
PNS/Torso Wearable Care Instructions	8-9

Post-Permanent FAQs	10-15
----------------------------------	-------



Quick Reference User Guide

Transmitter Assembly

Transmitter



Device is ON



Charging Error*

*A solid green light indicates a charging error. Turn the device OFF and beginning charging again.



Battery is Charging



Battery is Fully Charged





Antenna Disconnected



Device is OFF (no light)



Get Started with your device

1. Once the transmitter is placed in your wearable, place the antenna over the stimulator.
2. Press the transmitter ON button.  The green light will flash.
3. Press the  button ___ times.
4. Adjust position of antenna for best alignment and comfort.

- A** Twist Antenna Cable Here
- B** Increase
- C** Power On/Off
- D** Decrease
- E** Connect Charging Cable Here

Antenna

A solid yellow light indicates there is an issue between the antenna and the transmitter. Turn the device off and disconnect and reconnect the antenna to the transmitter. Turn the device back on. The yellow error light should now be gone and a flashing green light means the antenna and transmitter are communicating properly.

Charging

The transmitter **MUST** be turned off before charging. When plugged into an outlet, a steady blue light indicates charging. A flashing blue light indicates the transmitter is fully charged.

Press the **+** and the **-** buttons at the same time and release to change programs.

Program 1: Green light blinks ONE time

Program 2: Green light blinks TWO times

Program 3: Green light blinks THREE times

Quick Reference User Guide*

PNS Extremity Wearable



Quick Reference User Guide*

Torso Wearable



* Actual product placement will vary. Work with your physician and Curonix representative to identify proper placement. Product is to be worn over a thin layer of clothing at all times. Do not place directly on skin.



Wearable Care Instructions

For PNS Extremity and Torso

PNS Wearable Material:

75% Pa (Polyamide) / 20% Ea (Elastomer) / 5% PU (Polyurethane); Thread is 100% polyester; Velcro® is made by Velcro brand, and is an HTH model. The Curonix PNS Wearable Accessories (co-designed with Viviana Straps) use certified OEKO_TEX hypoallergenic fabric and gripper and are manufactured by hand in Italy.

PNS Wearable Care Instructions:

CAUTION: Remove all transmitter assembly components before washing. Machine wash at cold setting. Hang dry or tumble dry at a low temperature.

Torso Wearable Material:

Nylon, Polyacrylate (UBL) Polyethylene (hook), Polyester (Interfacing), Velcro® Brand Hook 805.

Torso Wearable Care Instructions:

CAUTION: Remove all transmitter assembly components before washing. Hand wash only with cold water; Wash with mild detergent; Air dry; Do not use bleach or other chemicals.

Transmitter Cover Material:

Polycarbonate/ABS blend.

Transmitter Cover Care Instructions:

Use water, isopropyl alcohol, or a mild solution of warm water and dish detergent and a soft cloth to clean the surface. Do not use harsh solvents or ammonia-based products to clean the cover.



For more information,
Please scan the QR Code or
visit us at our website.





Post - Permanent

FAQs

1. Will I have post-operative restrictions?

Your physician will provide discharge instructions. Common restrictions may include:

- Do not engage in rigorous physical activity such as twisting, bending, lifting heavy objects, or climbing.
- Do not remove or alter wound dressings without consulting clinical personnel.
- Continue to take any medications prescribed, such as antibiotics, unless otherwise directed by your clinician.
- Do not touch the wound or push, pull, or twist the implanted neurostimulator.

Temporarily, you may experience some pain around the implant site as the incisions heals. If you experience significant pain, swelling, discharge, or excessive redness around the wound, contact your clinician. In rare cases, adverse tissue reaction to implanted materials can occur during this period.

2. Will I receive a Curonix Medical Implant Card?

Yes, a Curonix Medical Implant Card will be mailed via FedEx to the address you provided within a month of the permanent procedure. The Curonix Medical Implant Card will contain information about your implanted device.

3. Who do I contact if I lose my Curonix Medical Implant Card?

Please contact Curonix Support at **(800) 965-5134** for a replacement.

Post - Permanent

FAQs

4. TSA / Airport Security

If you are wearing your Transmitter Assembly, you will need to remove it and turn it off prior to going through security. Place the Wearable and Transmitter Assembly on the security belt for scanning. Notify TSA that you have a medical implant and have your Curonix Medical Implant Card available if needed. They will direct you with any further instructions.

5. If I need an MRI, what do I need to do?

You will need to bring your Curonix Medical Implant Card with you and notify the facility that you have a medical device implanted. The Curonix Medical Implant Card contains the information needed for the facility to determine if they can perform the MRI.

6. What do I do if I move out of state and need programming or other device related support?

The local Curonix Representative can connect you to Curonix Support in your new area or you can call Curonix Support at **(800) 965-5134**.

Post - Permanent

FAQs

7. What do I do if I have to undergo a different medical procedure?

You should notify your health care provider(s) that you have a medical implant prior to receiving another medical procedure. Please refer your healthcare provider(s) to the Instructions for Use (IFU) available on curonix.com. The IFU includes specific information related to CT Scans, MRIs, radiofrequency ablations and other medical procedures.

8. Can I turn my Transmitter Assembly on while it is in the Wearable then place the Wearable on my body?

You may power on the Transmitter before placing it on your body, however, the Transmitter Assembly must be over the internal stimulator before increasing the Power Index.

9. How do I care for my Wearable?

PNS Extremity Wearable: Remove the Transmitter Assembly from the Wearable before washing on COLD. Hang dry or tumble dry on LOW temperature. Avoid heat as this will damage the wearable.

Torso Wearable: Remove the Transmitter Assembly from the Wearable before washing on COLD. **Air Dry**, avoid any heat as this will damage the Wearable.

You can access Wearable care instructions in the Transmitter Assembly Instructions for Use (IFU) available on curonix.com.

Post - Permanent

FAQs

10. How long does it take to recharge the Transmitter?

The Transmitter will typically recharge in 4-6 hours. If battery life is drained to depletion, recharging time may increase.

11. How often should I charge my Transmitter?

The Transmitter battery life varies depending on your individual therapy settings. We recommend charging the Transmitter every 12 hours if you are using the Transmitter 24/7.

12. Do I wear my Transmitter for the entire day?

Every patient is different; some require therapy on a continuous basis, while others will be able to wear it after a period of time for effective treatment. You will work with your physician to develop a therapy plan and the Curonix Representative will support your provider/physician's therapy plan.

13. Do I need to push the “plus” button on the Transmitter to my Power Index every time I put my Transmitter on?

Yes, the Transmitter automatically goes back to zero when turned off.

14. What if I miscount when pressing the “plus” button to my number?

You can turn the Transmitter off, then back on and push the plus button to your number.



Post - Permanent

FAQs

15. How long do I have to wait until I return to daily activities after the implant surgery?

Your physician will provide discharge instructions detailing your restrictions.

16. How long do I have to wait to shower?

Your physician will provide discharge instructions detailing your restrictions.

17. Can I get the Transmitter wet?

No, keep the Transmitter dry at all times.

18. My Transmitter won't turn on, what do I do?

Verify whether the Transmitter is turned off (no lights) before placing it on the charger. Allow the Transmitter to fully charge. If the issue persists, please contact your local Curonix Representative.

19. How do I know what program I am on when I turn the Transmitter on?

The center button will blink, indicating the Transmitter is on program 1, 2 or 3.

Post - Permanent

FAQs

19. How do I know what program I am on when I turn the Transmitter on? *(Continued)*



ON – Programming Setting 1
Blinking Green (Single)



ON – Programming Setting 2
Blinking Green (Double)



ON – Programming Setting 3
Blinking Green (Triple)

20. Can I go to sleep with my stimulator on?

Yes, you can sleep with the Transmitter powered on if needed.

21. Do I have to go through any follow-up appointments long-term with this device?

Check with your health care provider to schedule medical follow ups that may be needed. All medical questions will be answered by your health care provider. Contact your local Curonix Representative if you require any assistance or support in connection with the device.

22. Will the stimulator affect any of the medications I am currently on?

Please consult with your prescribing health care provider for any medication questions.

All medical questions should be directed to your health care provider.

CURONIX

The intent of this booklet is a reference guide, for more information on risks, warnings, and possible adverse side effects please refer to the Instructions for Use provided with the device or available on curonix.com.

All medical questions should be directed to your health care provider.

Curonix LLC is a manufacturer of medical devices and does not practice medicine. Only a physician can determine what treatment is appropriate. The contents of this document do not constitute medical, legal, or any other type of professional advice. For more information on risks, warnings, and possible adverse side effects please speak with your physician directly. For further product safety information refer to curonix.com.

For informational purposes only.